POLAR Horse Trainer

TRANSMITTER
Operating instruction
Gebrauchsanleitung
Mode d'emploi
Manual de instrucciones

POLAR®
heart rate monitors


English

Polar HorseTrainer Transmitter

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1. Moisten the horse's coat at the electrodes' position.
2. Identify the electrodes polarity:
The positive electrode **A** is marked (+) and is on the Polar logo side of the transmitter. The negative electrode **B** is marked (-).
3. Place the positive electrode on the left of the withers under the saddle or harness protection. Ensure that the grooves of the electrode is flat against the horse's skin. In some cases you may also attach the electrode around the girth by locking the clip.
4. Attach the plastic box **C** to the girth using the strap and clip. (For rider, attach the plastic box to the saddle cover)
5. Attach the negative electrode in a ventral position using the strap and clip.

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Español

Polar Horse Trainer Transmitter

HOW TO ATTACH THE POLAR HORSETRAINER TRANSMITTER

(See pictures for the correct location of the electrodes.)

1. Moisten the horse's coat at the electrodes' position.
2. Identify the electrodes polarity: The positive electrode **A** is marked (+) and is on the Polar logo side of the transmitter. The negative electrode **B** is marked (-).
3. Place the positive electrode on the left of the withers under the saddle or harness protection. Ensure that the grooves of the electrode is flat against the horse's skin. In some cases you may also attach the electrode around the girth by locking the clip.
4. Attach the plastic box **C** to the girth using the strap and clip. (For rider, attach the plastic box to the saddle cover)
5. Attach the negative electrode in a ventral position using the strap and clip.

Frequently asked questions

WHAT SHOULD I DO IF ...

... THE HEART RATE SEEMS TO BE WRONG?

- Make sure that the positions of the electrodes are correct.
- Make sure that the horse's coat under electrodes is thoroughly moistened.
- Accumulated sweat and dirt may interfere the measurement of the signal from the heart to the transmitter. Check if the electrodes need cleaning.
- A lot of using could cause oxidization of the electrodes. Wash the electrodes with water and soap solution, then dry them carefully.
- Make sure that the cables do not move too much.
- In some cases (at about 120 bpm) a double pulse could occur at each heart beat (R + T waves). This problem disappears when the intensity increases.
- If you are in an electromagnetic field (e.g. near high voltage power lines), move further away.

... THERE ARE TRANSMISSION CUTS?

- Make sure that the transmitter is in a right position.
- In rider application, check that the watch receiver is inside the transmission range and not further than 1metre/3 feet from the transmitter.
- The transmission distance may have decreased. Have the battery checked by authorised Baumann & Haldi SA service personnel. Contact your distributor/dealer for details.

... THE WATCH RECEIVER DISPLAY IS BLANK OR FADING?

Have the battery checked by authorised Baumann & Haldi SA service personnel. Contact your distributor/dealer for details.

Care and maintenance

TRANSMITTER

The Polar HorseTrainer transmitter is active only when the electrodes are in contact with the horse's coat. Sweat and moisture can keep the electrodes wet and the transmitter activated, which shortens the battery life. Therefore, clean the electrodes and wipe them dry after use.

Make sure that the electrode (conductive) are not in contact.

Dirt on the electrodes weakens the contact between a horse's body and the electrodes. To ensure good operation, wash the electrodes regularly with a mild soap and water solution. However, never use alcohol or any abrasives such as steel wool or clening chemicals. They will cause permanent damage to the electrodes.

The material used for electrodes and strap is very soft and comfortable, take care of these parts to prevent any damage. Electrodes should be replaced by an authorised Baumann & Haldi SA service personnel. For details, contact your dealer or distributor.

BATTERIES

The transmitter's estimated average battery life is 1000 hours.

Battery replacement is recommended to be done by service personnel authorised by Baumann & Haldi SA. For details, contact your dealer or distributor.